

## Get the right fit with the CervAlign®.

- CervAlign's patented pillow fitting system helps to ensure the best fit possible.
- Provides comfort and support through a full range of treatment.
- Curved cervical lobe follows the neckline, maintaining proper cervical support in normal sleeping positions.
- Doubles as a conventional pillow for normal sleeping comfort or to help your body adjust to sleeping on the cervical lobe.
- For back and side sleepers.
- Patent #6,471,661



## Support your neck in its natural position.

CervAlign's unique design and custom fit provide proper cervical support while sleeping on your back or side. A proper sleeping posture can provide relief from, and help prevent:

- Tension headaches
- Neck muscle and joint sprains and strains
- Ligament sprains
- Whiplash injuries
- Arthritis
- Tingling and numbness in the arms and hands

To order, contact:



Manufactured by:  
**Core Products**  
INTERNATIONAL, INC.  
800-365-3047 • 715-294-2050 (Int'l)  
[www.coreproducts.com](http://www.coreproducts.com)

# People come in different sizes.



**Core Products**  
INTERNATIONAL, INC.  
[www.coreproducts.com](http://www.coreproducts.com)

*Making Your Life More Comfortable*

# Experience the CervAlign® difference.

**NEW FEATURE!**  
Sewn-in Indented Center  
provides better neck  
support for back sleeping.

Large upper pillow section  
provides normal sleeping  
comfort when the neck is  
not stiff or irritated.

Made with 230 thread-count fabric and new  
generation fiber fill for maximum comfort.



Curved cervical lobe follows  
the neckline, maintaining proper  
neck support as you roll from  
your back to your side.

Curved design allows for  
shoulder comfort with mild  
neck extension during sleep.



Correct Fit



Incorrect Fit



## Start with a proper fit.

Choosing a correctly sized pillow is a vital step in proper neck care. An improperly sized corrective pillow will not provide the proper support to return your spine to a healthy position, which in many cases will prolong the period of treatment. The CervAlign Pillow Fitting Guide helps to ensure the most rapid recovery possible.

### Fitting guides allow for selection of proper pillow size:

1. Have the individual stand up straight.
2. Place the rounded edge of the sizing guide to the nape of the neck.
3. Visualize a vertical line extending upward from the edge of the shoulder blades. The flat edge of the sizing guide should land on this line.

If you are having trouble visualizing the line, let the remaining guides hang down against the individual's back. If the guides rest gently against the back, you have found the correct fit. If the guides hang out away from the body, try the next size down.

## Unique design and fit.

The CervAlign is different from most pillows – it is two pillows in one! It features a unique cervical lobe, fitted just for you. This curved lobe places your neck in its most natural, neutral position and maintains this proper cervical support as you roll from your back to your side. It is the perfect complement to your professional treatments. Also, when your neck is not stiff or irritated, you can turn the pillow around and use the large upper pillow section for normal sleeping comfort.

The CervAlign is made with 230 thread-count fabric and new generation fiber fill for maximum comfort. Washable. Available in three sizes: 5", 6" and 7". See your healthcare professional to determine the fit that's right for you.

Pillow fitting guides allow for  
selection of proper lobe size.  
Patent# 6,471,661



Product #265, 5" Lobe  
Product #266, 6" Lobe  
Product #267, 7" Lobe